

TRADITIONAL *MARLEN Ostomy Products*

GENERAL MEDICAL PRODUCTS

CONTINENT OSTOMY SUPPLIES


A complete range of catheters and accessories for continent ostomies.

Curved Catheters

Large No. 34 French
Total length — 18" (46 cm)



Medium No. 30 French
Total length — 18" (46 cm)



Straight Catheters


Large No. 34 French
Total length — 15" (38 cm)



Medium No. 30 French
Total length — 15" (38 cm)



Small No. 24 French
Total length — 15" (38 cm)



Catheter Support


For large catheter — 15070

For medium catheter — 15080



Elastic Belt


(For catheter support)



15090

Leg Bag Kit (Complete)

15100



Catheter Connectors


Large — 15140



Medium/Small — 15150



15110



Bedside Collector Kit (Complete)

Syringe (60 cc)



15120



Marlen Manufacturing & Development Co.
5150 Richmond Road • Bedford, Ohio 44146 • USA
Phone (216) 292-7060 • Fax (216) 292-9196
E-mail: info@marlenmfg.com • Web: www.marlenmfg.com

TRADITIONAL *MARLEN Ostomy Products*

GENERAL MEDICAL PRODUCTS

Loop Ostomy Rods

Intended for assistance in maturing loop ostomies. Made of durable plastic. Disposable one-time use.



616 — 2" long with eyelets at both ends



617 — 3" long with eyelets at both ends



618 — 3" long with eyelet at one end and T-Bar at other end



619 — 4" long with eyelets at both ends




620 — 4" long with eyelet at one end and T-Bar at other end

Biliary Drainage Kit


BK-6

Ideal for biliary and other drainages.

Lightweight pouch features a "T-Tube" connector and built-in flutter valve to prevent backflow. Push-Pull valve for draining. One size. Elastic waist belt also included.



Available in clear.



Length "A"	Width "B"	Capacity
9-1/4 in./23.5 cm	5-1/4 in./13.3 cm	20 oz./590 cc

Colonic Enema Set

MCE-850

For safe, effective evacuation in cases of decreased sphincter control. Also ideal for barium enemas. Irrigating bag holds 48 oz./1500 cc.

