


TRADITIONAL *MARLEN Ostomy Products*

ACCESSORIES



Part No. 9000

MARLEN Skin Shield™


Protective Adhesive Skin Barrier is non-absorbent and forms a long-lasting seal. Package of five. Each pad is 4" x 4" (10cm x 10cm) with 1/2" (12mm) starter opening. May be ordered without starter hole.



Part No. P-116

MARLEN Protex™ Powder (Gum Karaya)


Used to protect the skin against excoriation from ostomy drainage. 4 oz. (120cc) bottle.



Part No. 9100

MARLEN AquaTack™


Protective Hydrocolloid Skin Barrier adheres on dry or moist skin. Resists drainages. Package of five. Each pad is 4" x 4" (10cm x 10cm) with 1/2" (12mm) starter opening. May be ordered without starter hole.



Part No. PS-176

MARLEN Protex™ Sheet


Ideal protection for irregular shaped stomas and hard-to-fit areas. Individually packaged. 8" (20.3cm) square. Made of gum karaya. Custom shapes and sizes can be easily cut.



Part No. 107
Part No. 107-0

MARLEN Double-Faced Adhesive Tape Discs

Used to adhere traditional pouch/mounting ring appliances to body. Package of 10. 3-7/8" outside diameter (9.9cm). Pre-cut stoma openings up to 1-5/8" (4.2cm).



Part No. FP-601

MARLEN Plastic Foam Pads

Used to help build up the areas around surgical scars and depressions to provide better adherence. Package of 10. 4" (10.2cm) diameter with pre-cut stoma openings.



Part No. XLF-171

MARLEN Extra-Large Double-Faced Adhesive Tape Discs


Shaped to accommodate the Extra-Large XTL-55 and XTL-75 mounting rings. The outside dimensions are 3-1/2" x 4-1/2" (8.9cm x 11.5cm). Pre-cut stoma openings of 1-3/4" to 2-1/2" (4.4cm to 6.4cm).



Part No. 103


MARLEN Rubber O-Ring Seal

Individually packaged. Used with all mounting rings (except Extra-Large XTL Series) to assemble the mounting ring to the pouch.



Part No. XTL-161

Individually packaged. Used only with Extra-Large XTL-55 and XTL-75 Mounting Rings.



Part No. 112

MARLEN Small Protex™ Powder Pad (Karaya Seal)


Worn under appliance for protection from drainage coming in contact with the skin. Package of 10. Pre-cut with stoma opening sizes up to 1-1/2" (3.8cm).



Part No. 605

MARLEN Lock Ring

Individually packaged. Can be used as an alternative to Rubber O-Ring Seal with All Flexible Mounting Rings.



Part No. PP-C-175

MARLEN Large Protex™ Powder Pad (Karaya Seal)

Package of 5. Outside diameter of 3-3/4" (9.5cm) will accommodate openings up to 2-1/2" (6.4cm).




Part No. 616

MARLEN Rubberbands

Package of 25. Used to attach Urinary Drain Fitting (Pt. No. 108) to pouch.

TRADITIONAL *MARLEN Ostomy Products*

ACCESSORIES



Part No. 104

MARLEN Elastic Appliance Belt (Adult)

Used with traditional mounting rings for additional support. Individually packaged.

Part No. 105

MARLEN Elastic Appliance Belt (Child)


Individually packaged.



Part No. 5004

MARLEN Ultra™ Elastic Waist Belt

Used with Ultra™ one-piece and Duet two-piece appliances for additional support. Individually packaged.



Part No. 117

MARLEN Plastic Clamp Closure


Pouch closure for re-usable ileostomy (drainable) pouches. Individually packaged.



Part No. 118

MARLEN "Snap-Lok"™ Pouch Closure

One-piece pouch closure for re-usable ileostomy (drainable) pouches. After rolling up spout, sides fold and snap to secure. Individually packaged.



Part No. 5018

MARLEN Ultra™ "Snap-Lok"™ Pouch Closure

One-piece pouch closure designed for the wider, Ultra™ drain spout. After rolling up spout, sides fold and snap to secure. Individually packaged.



Part No. 5019

MARLEN Ultra™ "Comfort Clip"™ Pouch Closure

Scissor-like action securely seals bottom of wider Ultra™ drain spout. Individually packaged.



Part No. 108

MARLEN Urinary Drainage Fitting

Used as drain fitting for MAF-12 and MU-3R pouches. Attach with MARLEN Rubberband (Pt. No. 606). Individually packaged.



Part No. 310

MARLEN Screw Connector

Used to attach MAF-12 and MU-3R pouches to night drain bags or bottles and leg bags. Individually packaged.



Part No. 603

MARLEN Stoma Guide Strips

Dissolving paper strips used to help guide appliance over stoma. Package of 100.



Part No. 503

MARLEN Liquid "Odour-Guard"™

4 oz. (120cc) bottle
Used to neutralize odors in the pouch.



Part No. 504

MARLEN Deodorizer

12 oz. (360cc) bottle
Ideal for cleaning re-usable pouches and mounting rings. Concentrated — makes up to 6 gallons of solution.



Part No. 507

MUCOSPERSE™ Mucous Dispersant

4 oz. (120cc) bottle
Dissolves mucous build-up in urinary pouches and collection devices.



Part No. 402

MARLEN Plastic Squeeze Bottle


1 oz. (30cc) bottle. Handy travel size for dispensing Odour-Guard, Deodorizer, and Mucosperse.



Part No. 106

MARLEN Plastic Airing Hook

Convenient for hanging appliance to dry after cleaning. Individually packaged.



Part No. 612

3M Micropore Tape

Lightweight "paper" tape (1" width) used for additional support. Package of 3 rolls.



Marlen Manufacturing & Development Co.

5150 Richmond Road • Bedford, Ohio 44146 • USA

Phone (216) 292-7060 • Fax (216) 292-9196

E-mail: info@marlenmfg.com • Web: www.marlenmfg.com