

TRADITIONAL *MARLEN Ostomy Products*

"ALL FLEXIBLE" MOUNTING RINGS

- Accommodates stoma openings of 1/2" to 1-5/8" (1.27cm to 4.13cm) in 1/16" increments.

- Compounded of green neoprene or white vinyl
- Soft, flexible and completely pliable

FLAT "All-Flexible" Mounting Rings

- For persons with flat abdomens, good muscle tone and protruding stomas




Basic All-Flexible

- For average or larger build with centrally-located stoma and adequate clearance on all sides. 3-3/4" (9.5cm) diameter

Part No. WV-101 White Vinyl

Part No. GN-101 Green Neoprene Rubber




Oval Flat

- For stomas located closer to the waistline where bending may be an important factor. 3-3/4" x 2-3/4" (9.5cm x 7.0cm)

Part No. WV-50 White Vinyl

Part No. GN-50 Green Neoprene Rubber



Compact Flat

- For slimmer persons with smaller torsos. 3-3/4" x 3-1/8" (9.5cm x 7.9cm)

Part No. WV-45 White Vinyl

Part No. GN-45 Green Neoprene Rubber



Petite Flat


- For the small person with a flat abdomen or a stoma in a hard-to-fit area close to the rib cage or hips. 3-3/4" x 2-7/8" (9.5cm x 7.3cm)

Part No. WV-40 White Vinyl

Part No. GN-40 Green Neoprene Rubber

CONVEX "All-Flexible" Mounting Rings

- For persons with fleshier abdomens and short or flush stomas. The 1/4" (6mm) convexity applies slight pressure to the stomal area, allowing the stoma to stay within the mounting ring opening and not retract.



Oval Convex

- For stomas located close to the waistline where bending may be an important factor. 3-3/4" x 2-3/4" (9.5cm x 7.0cm)

Part No. WV-60 White Vinyl

Part No. GN-60 Green Neoprene Rubber



Compact Convex

- For the smaller person with a well-positioned stoma. 3-3/4" x 3-1/8" (9.5cm x 7.9cm)

Part No. WV-102 White Vinyl

Part No. GN-102 Green Neoprene Rubber



Petite Convex

- For the small person with soft abdomen or with a stoma in a hard-to-fit area, close to the rib cage or hips. 3-3/4" x 2-7/8" (9.5cm x 7.3cm)


Part No. WV-42 White Vinyl

Part No. GN-42 Green Neoprene Rubber

DEEP CONVEX "All-Flexible" Mounting Rings

- Soft and flexible with the deepest convexity in a mounting ring

- Accommodates stoma openings of 1/2" to 1-5/8" (1.27cm to 4.13cm)
- Compounded of green neoprene or white vinyl



Basic Deep Convex

- For flabby abdomens and short or flush stomas where deep convexity is needed and an all-flexible ring is desired. Deepest convexity in a mounting ring — 3/8" (1cm). 3-3/4" (9.5cm) diameter

Part No. WV-80 White Vinyl

Part No. GN-80 Green Neoprene Rubber



Compact Deep Convex

- All-flexible comfort with reduced diameter to fit short or flush stomas positioned close to rib cage, hip or navel or for the shorter torso. Deepest convexity in a mounting ring — 3/8" (1cm). 3-3/4" x 3-1/8" (9.5cm x 7.9cm)

Part No. WV-85 White Vinyl


Part No. GN-85 Green Neoprene Rubber

TRADITIONAL *MARLEN Ostomy Products*

"ALL-FLEXIBLE" MOUNTING RINGS

EXTRA LARGE "All-Flexible" Mounting Rings

- Soft, flexible and completely pliable
- Designed specifically to accommodate larger stoma openings and double-barreled stomas of 1-3/4" to 2-1/2" (4.45cm to 6.35cm)
- Oval shape allows maximum adhesive area without added bulk
- Compounded of white vinyl



"XTL" Oval Flat

- For persons with flat abdomens and good muscle tone. 4-1/2" x 3-1/2" (11.4cm x 8.9cm)

Part No. XTL-55 White Vinyl




"XTL" Oval Convex

- Firmer support for the person with a fleshier abdomen, lacking muscle tone. Deep convexity—3/8" (1cm). 4-1/2" x 3-1/2" (11.4cm x 8.9cm)

Part No. XTL-75 White Vinyl

"SEMI-FLEXIBLE" MOUNTING RINGS

- Firmer mounting ring designed for added support in the stomal area
- Constructed with a stainless steel core covered with soft neoprene rubber
- Flexible outer edges for comfort
- Accommodates stoma openings of 1/2" to 1-5/8" (1.27cm to 4.13cm) in 1/16" increments.



Basic Regular Convex

- Ideal for slight or medium build persons, especially in the immediate post-operative period. Very slight convexity—7/32" (5mm). 3-5/8" (9.2cm) diameter

Part No. SF-10



Basic Medium Convex

- For stouter persons or for short, flush stomas, this ring offers maximum adhesive area for stomas positioned to accommodate it. Medium convexity—9/32" (7mm). 4" (10.2cm) diameter

Part No. SF-15



Compact Regular Convex

- Smaller in circumference with convexity of 7/32" (5mm) to accommodate short stomas. 3-5/8" x 3-1/8" (9.2cm x 7.9cm)

Part No. SF-20



Compact Medium Convex

- Reduced in diameter to fit short or flush stomas positioned close to rib cage, hip or navel or for the shorter torso. Medium convexity—9/32" (7mm). 3-5/8" x 3-1/8" (9.2cm x 7.9cm)

Part No. SF-65



MARLEN

Marlen Manufacturing & Development Co.

5150 Richmond Road • Bedford, Ohio 44146 • USA

Phone (216) 292-7060 • Fax (216) 292-9196

E-mail: info@marlenmfg.com • Web: www.marlenmfg.com